



"Improving Patient Care Through Research & Education"

Home > Clinical Evaluations > Evaluations

ESTELITE FLOW QUICK (2010 PPL) ++++1/2

EDITORS' EVALUATION

Tokuyama Dental Corp./Tokuyama America Inc.
877.378.3548
www.tokuyama-us.com

Editors' Comments

- ▶ "Good esthetics and finishing characteristics."
- ▶ "Controlled flow and a lack of bubbles."
- ▶ "Increase the opacity of the opalescent shades."

EXCLUSIVE REVIEW



Featured in:
January-February
2010 Volume 27, No.
01

SHARE

Subscriber Toolbox

[READ FULL ISSUE](#)
[ADD TO MY ARCHIVE](#)
[PRINT THIS PAGE](#)

More Information

[Winner of the 2010 Top Flowable Composite Award](#)

Description

ESTELITE FLOW QUICK is a low-viscosity, light-cured, radiopaque, submicron-filled resin composite. It contains 71% by weight (53% by volume) silica-zirconia and silica-titania filler. All filler particles are spherical (mean particle size 0.3 μm , particle range 0.04-0.6 μm), designed to provide high retention of gloss and high wear resistance. **ESTELITE FLOW QUICK** is indicated for direct placement in anterior and posterior restorations, block out of undercuts, and repair of ceramic/composite restorations. The prominent feature of **ESTELITE FLOW QUICK** is RAP (Radical Amplified Photopolymerization) technology that promotes more rapid polymerization. As reported by the manufacturer, curing time can be shortened by one-third with curing lights having a wavelength range of 400-500 nm (peak 470 nm). Each consultant received one, 1-ml syringe each of A1, A2, A3, A3.5, B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, BW, incisal, and cervical of **ESTELITE FLOW QUICK**; syringe tips; and instructions. Fifteen editors evaluated this composite in over 460 restorations. This product received a 93% clinical rating.

Suggested Retail Cost

- \$26.00/syringe
- \$46.80/2-syringe pack
- \$70.00/3-syringe pack

Evaluation Protocol

- Fifteen editors of THE DENTAL ADVISOR clinically evaluated **ESTELITE FLOW QUICK** over a six-week period. 426 restorations were placed.
- **ESTELITE FLOW QUICK** was evaluated on a 1-5 rating scale:
1 = poor, 2 = fair, 3 = good, 4 = very good, 5 = excellent.

Editors' Observations

Manufacturer's Instructions/Packaging

Instructions are clearly printed with the curing times for the variety of shades and the types of curing lights frequently used. The syringe is not bulky and is easily manipulated. Each syringe is clearly labeled for quick shade recognition, but the editors noted that more material in each syringe would be an improvement.

Ease of Dispensing/Lack of Bubbles

Syringes come with plastic application tips that are curved to allow easy placement in a variety of locations. The material is easily expressed from the syringe. The editors found little or no presence of bubbles.

Viscosity/Stackability/Wettability

ESTELITE FLOW QUICK has a higher viscosity than most flowable materials. It flows well into the preparation, yet has good stackability. There is a sacrifice in wettability with this material due to the high viscosity and stackability. The editors reported they would select *PALFIQUE ESTELITE LV High Flow*, a lower-viscosity material, in areas where wettability was desired.

Shade Match to Vita Shade Guide

The editors found the shades matched the Vita Shade Guide adequately. They noted that 17 available shades were too many; adding an opaquer for blocking is suggested.

Blending/Translucency/Opacity

ESTELITE FLOW QUICK blends well due to excellent translucency, allowing for a good chameleon effect and easy addition in a layering technique for anterior restorations. An increase in the opacity of the opalescent shades is recommended.

Radiopacity

Radiographic appearance of **ESTELITE FLOW QUICK** is very close to that of dentin.

Sensitivity to Ambient Light

There was no report of sensitivity to ambient light. The RAP technology did not impair workability and was reported as one of the most liked features by the editors.

Finishing and Polishing

The spherical particle shape and small size allow for easy finishing and polishing. The editors reported that the final polish was very good.

Final Restoration

Overall, the editors were satisfied with the qualities of **ESTELITE FLOW QUICK**. Though not a typical flowable composite, the material has very good handling properties, blends well, cures easily, and has an excellent overall finish that makes a quality restoration.

Clinical Tip

- To encourage good adhesion, the manufacturer recommends not using a dual- or self-curing bonding agent with **ESTELITE FLOW QUICK**.

More About: [Tokuyama Dental Corp./Tokuyama America Inc.](#),
[Composite - Flowable, Composite, Flowable Composite](#)

Editors' Note: *The original evaluation was published in THE DENTAL ADVISOR, Vol. 23, No. 3, April 2006.*

CE0981

[FAQ](#) [Copyright Policy](#) [Privacy Policy](#) [Terms & Conditions](#) [Contact](#)

Dental Consultants, Inc. 3110 West Liberty, Ann Arbor, Michigan 48103
(800) 347-1330 (734) 665-2020 (734) 665-1648 fax info@dentaladvisor.com
© 2009 Dental Consultants, Inc.